

**A rapid test for the qualitative detection of Synthetic Marijuana in human urine.  
For medical and other professional in vitro diagnostic use only.**

**INTENDED USE**

The **DrugControl K2 Test** is a rapid chromatographic immunoassay for the detection of Synthetic Marijuana metabolite in human urine. The synthetic marijuana detected by the test includes, but are not limited to, the metabolites of JWH-018 and JWH-073. The following table lists compounds that are positively detected in urine by the **DrugControl K2 Test** at 5 minutes:

TEST DEVICE	CALIBRATOR / related compounds	CUT-OFF LIMIT VALUE [ng / ml]
K2 50	JWH-073 4-butanoic acid metabolite	50
	JWH-018 5-Pentanoic acid metabolite	50
	JWH-018 4-Hydroxypentyl metabolite	400
	JWH-018 5-Hydroxypentyl metabolite	500
	JWH-073 4-Hydroxybutyl metabolite	500

This assay provides only a qualitative, preliminary, analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

**SUMMARY**

Synthetic Marijuana or K2 is a psychoactive herbal and chemical product that, when consumed, mimics the effects of Marijuana. It is best known by the brand names K2 and Spice, both of which have largely become genericized trademarks used to refer to any synthetic Marijuana product. The studies suggest that synthetic marijuana intoxication is associated with acute psychosis, worsening of previously stable psychotic disorders, and also may have the ability to trigger a chronic (long-term) psychotic disorder among vulnerable individuals such as those with a family history of mental illness.

Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 72 hours after smoking (depending on usage/dosage).

As of March 1, 2011, five cannabinoids, JWH-018, JWH-073, CP-47, JWH-200 and cannabicyclo hexanol are now illegal in the US because these substances have the potential to be extremely harmful and, therefore, pose an imminent hazard to the public safety.

The **DrugControl K2 Test** is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of synthetic marijuana metabolite in human urine. The **DrugControl K2 Test** yields a positive result when the synthetic marijuana metabolite in urine exceeds 50 ng/mL.

**TEST PRINCIPLE**

The **DrugControl K2 Test** is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Synthetic Marijuana metabolite, if present in the urine specimen below 50ng/mL, will not saturate the binding sites of antibody-coated particles in the test device. The antibody-coated particles will then be captured by immobilized synthetic marijuana conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Synthetic Marijuana metabolite level exceeds 50ng/mL because it will saturate all the binding sites of anti- Synthetic Marijuana antibodies.

A drug-positive urine specimen will not generate a colored line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

**REAGENTS**

The test strip contains mouse monoclonal anti- synthetic marijuana antibody-coupled particles and synthetic marijuana-protein conjugate. A goat antibody is employed in the control line system.

**PRECAUTIONS**

- For medical and other professional in vitro diagnostic use only.
- Do not use after the expiration date.
- The test device should remain in the sealed pouch until use.
- Do not use the test if the foil pouch is damaged
- Do not moisten nitrocellulose membrane with urine samples.
- Read the entire procedure carefully prior testing.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- Humidity and temperature can adversely affect results.
- The used test device should be discarded according to federal state and local regulations.
- Do not reuse tests
- Avoid cross-contamination of urine samples by using a new specimen collection container for each urine sample.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.

**STORAGE AND STABILITY**

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Tests should be kept out of direct sunlight.

- Do not freeze
- Do not use beyond the expiration date.

**SPECIMEN COLLECTION AND PREPARATION**

**Urine Assay**

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain clear specimen for testing.

**Specimen Storage**

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

**MATERIALS PROVIDED**

- Test in pouch
- Disposable sample dropper (in pouch)
- Package insert


**MATERIALS REQUIRED, BUT NOT PROVIDED**

- Specimen collection container
- Timer
- Positive and negative controls

**DIRECTIONS FOR USE**

1. Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.
2. Bring the pouch to room temperature before opening it.
3. Remove the Test from the sealed pouch and use it within one hour.
4. Place the test cassette on a clean and level surface.
5. Hold the dropper vertically and transfer 3 full drops of urine to the specimen well (S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well (S).
6. Read results at 5 minutes. Do not interpret the result after 10 minutes.

**INTERPRETATION OF RESULTS**



- Positive:** One color line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the K2 concentration is above the detectable cut-off level. (substances & cut-off concentrations see table on page 1).
- Negative:\*** Two lines appear. One color line should be in the control region (C), and another apparent color line should be in the test region (T). This negative result indicates that the K2 concentration is below the detectable cut-off level.
- Invalid:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test immediately and contact your local distributor.
- \* **Note:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

**QUALITY CONTROL**

A procedural control is included in the test. A color line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test; however it is recommended that positive and negative controls be tested as good laboratory testing practices to confirm the test procedure and to verify proper test performance.

**LIMITATIONS**

- The **DrugControl K2 Test** provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrophotometry (GC/MS) is the preferred confirmatory method<sup>1,2</sup>.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- This test does not distinguish between drugs of abuse and certain medications.
- The **DrugControl K2 Test** is intended for use with human urine specimens only.

**EXPECTED VALUES**

This negative result indicates that the synthetic marijuana metabolite concentration is below the detectable level of 50ng/ml. Positive result means the concentration of synthetic marijuana metabolite is above the level of 50ng/ml. The **DrugControl K2 Test** has a sensitivity of 50ng/ml

**PERFORMANCE CHARACTERISTICS**

**Accuracy**

A side-by-side comparison was conducted using the **DrugControl K2 Test** and GC/MS. The following results were tabulated:

<b>K2</b>		<b>GC/MS</b>		<b>Total Results</b>
		<i>Positive</i>	<i>Negative</i>	
<b>DrugControl K2 Test</b>	<i>Positive</i>	78	3	81
	<i>Negative</i>	2	167	169
<b>Total Results</b>		80	170	250
<b>% Agreement with this Test</b>		<b>97.5%</b>	<b>98.2%</b>	<b>98.0%</b>

**Analytical Sensitivity**

A drug-free urine pool was spiked with K2 at the following concentrations: 0 ng/mL, 25 ng/mL, 37.5 ng/mL, 50 ng/mL, 62.5 ng/mL, 75 ng/mL and 150 ng/mL. The results demonstrate >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

<b>K2 Concentration (ng/mL)</b>	<b>Percent of Cut-off</b>	<b>n</b>	<b>Visual Result</b>	
			<i>Negative</i>	<i>Positive</i>
0	0	30	30	0
25	-50%	30	30	0
37.5	-25%	30	26	4
50	<i>Cut-off</i>	30	15	15
62.5	+25%	30	3	27
75	+50%	30	0	30
150	3X	30	0	30

**Precision**

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to HPLC, no synthetic marijuana, 25% synthetic marijuana above and below the cut-off, and 50% synthetic marijuana above and below the 50ng/mL cut-off was provided to each site. The following results were tabulated:

<b>K2 Concentration (ng/mL)</b>	<b>n per Site</b>	<b>Site A</b>		<b>Site B</b>		<b>Site C</b>	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
25	10	10	0	10	0	10	0
37.5	10	8	2	8	2	9	1
62.5	10	1	9	2	8	2	8
75	10	0	10	0	10	0	10

**Effect of Urinary Specific Gravity**

Fifteen (15) urine samples of normal, high, and low specific gravity ranges were spiked with 25ng/mL and 75ng/mL of synthetic marijuana. The **DrugControl K2 Test** was tested in duplicate using the fifteen neat and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

**Effect of Urinary pH**

The pH of an aliquoted negative urine pool was adjusted to a pH range of 4 to 9 in 1 pH unit increments and spiked with synthetic marijuana to 25ng/mL and 75ng/mL. The spiked, pH-adjusted urine was tested with the **DrugControl K2 Test** in duplicate. The results demonstrated that varying ranges of pH do not interfere with the performance of the test.

**CROSS-REACTIVITY**

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or synthetic marijuana positive urine. The following compounds show no cross-reactivity when tested with the **DrugControl K2 Test** at a concentration of 100 µg/mL.

### Non Cross-Reacting Compounds

4-Acetaminophenol	R (-)Deprenyl	3-Hydroxytyramine	Nifedipine	Serotonin
Acetone	Dextromethorphan	(Dopamine)	Nimesulide	(5-Hydroxytryptamine)
Acetophenetidin	Diazepam	Hydroxyzine	Norcodeine	Sodium chloride
N-Acetylprocainamide	Diclofenac	Ibuprofen	Norethindrone	Sulfamethazine
Acetylsalicylic acid	Dicyclomine	Imipramine	d-Norpropoxyphene	Sulindac
Albumin	Diflunisal	Iproniazide	Noscapine	Sustiva (Efavirenz)
Amitriptyline	4-Acetaminophenol	(-)-Isoproterenol	d,l-Octopamine	Temazepam
Amobarbital	Acetone	Isoxsuprine	Orphenadrine	Tetracycline
Amoxapine	Acetophenetidin	Kanamycin	Oxalic acid	Tetrahydrocortexolone
Amoxicillin	4-Dimethylaminoantipyrine	Ketamine	Oxazepam	Tetrahydrocortisone,
Ampicillin	Diphenhydramine	Ketoprofen	Oxolinic acid	3-acetate
Ascorbic acid	5,5-Diphenylhydantoin	Labetalol	Oxycodone	Tetrahydrozoline
Aminopyrine	Disopyramide	Levorphanol	Oxymetazoline	Thebaine
Apomorphine	Doxylamine	Lidocaine	Oxymorphone	Thiamine
Aspartame	Ecgonine	Lindane	Papaverine	Thioridazine
Atropine	Ecgonine methylester	(Hexachlorocyclohexane)	Pemoline	l-Thyroxine
Benzilic acid	EMDP	Loperamide	Penicillin-G	Tolbutamide
Benzoic acid	Ephedrine	4-Dimethylaminoantipyrine	Pentazocine	cis-Tramadol
Benzphetamine	l-Ephedrine	Diphenhydramine	Perphenazine	trans-2-
Bilirubin	l-Epinephrine	5,5-Diphenylhydantoin	Phencyclidine	Phenylcyclopropylamine
Brompheniramine	(±)-Epinephrine	Maprotiline	Phenelzine	Trazodone
Buspirone	Erythromycin	Meperidine	Pheniramine	Trimethobenzamide
Cannabinol	β-Estradiol	Meprobamate	Phenobarbital	Triamterene
Cimetidine	Estrone-3-sulfate	d-Methamphetamine	Phenothiazine	Trifluoperazine
Chloral hydrate	Ethanol (Ethyl alcohol)	l-Methamphetamine	Phentermine	Trimethoprim
Chloramphenicol	Ethyl-p-aminobenzoate	Methadone	Prednisolone	Trimipramine
Chlordiazepoxide	Etodolac	Methoxyphenamine	Prednisone	Tryptamine
Chloroquine	Famprofazone	(+)-3,4-Methylenedioxy-	Maprotiline	d,l-Tryptophan
Chlorothiazide	Fentanyl	Methylphenidate	Meperidine	Tyramine
(+)-Chlorpheniramine	Fluoxetine	Mephentermine	Meprbamate	d,l-Tyrosine
(±)-Chlorpheniramine	Furosemide	Metoprolol	Procaine	Uric acid
Chlorpromazine	Genitisc acid	Morphine-3-β-D-	Promazine	Verapamil
Chlorprothixene	d-Glucose	glucuronide	Promethazine	Digoxin
Cholesterol	Guaiacol glyceryl ether	Morphine sulfate	l-Propoxyphene	Lithium carbonate
Clomipramine	Hemoglobin	Methyprylon	d,l-Propranolol	l-Phenylephrine
Codeine	Hydralazine	Nalidixic acid	d-Pseudoephedrine	Procaine
Cortisone	Hydrochlorothiazide	Nalorphine	Quinacrine	Promazine
(-)-Cotinine	Hydrocortisone	Naloxone	Quinidine	Promethazine
Creatinine	o-Hydroxyhippuric acid	Naltrexone	Quinine	
Cyclobarbitol	p-Hydroxymeth-	α-Naphthaleneacetic acid	Ranitidine	
Cyclobenzaprine	amphetamine	Naproxen	Riboflavin	
Deoxycorticosterone		Niacinamide	Salicylic acid	

### LIMITATIONS

It is impossible to check any and all - other than those drugs mentioned in the product insert - for cross-reactivity or any other influences to the to be detected drug of abuse ( DOA ).

If the patient takes a „cocktail“ of several different drugs or medication cannot be excluded that a non-reproducible cross-reaction can falsified the test result.

### BIBLIOGRAPHY

1. Wong, R., the Current Status of Drug Testing in the US Workforce, Am.Clin.Lab, 2002; 21(1):21-23
2. Info Facts -Club drugs, NIDA, May 2006, <http://www.nida.nih.gov/infofacts/clubdrugs.html>
3. Drug Fact Sheet, DEA, January 2012, <http://www.dea.gov>
4. Wong, R., the Effect of Adulterants on Urine Screen for Drugs of Abuse: Detection by an On-site Dipstick Device, Am.Clin.Lab, 2002; 21(3); 14-18
5. Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man, Biomedical Publications, Davis, CA, 1982.
6. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.

	Manufacturer		Contents sufficient for <n> tests
	For in vitro diagnostic use only		Lot. no.
	For single use only		Expiration date
	Read instructions for use		Store at
	Keep away from direct sunlight		Ordering number
	Keep dry		

This operating manual conforms to the latest technology / revision. Subject to change without prior notice!



ulti med Products (Deutschland) GmbH  
Reeshoop 1 • 22926 Ahrensburg • Germany  
Telefon: +49-4102 – 80090  
Fax: +49-4102 – 50082  
e-mail: [info@ultimed.de](mailto:info@ultimed.de)

Distributor:  
ulti med Products (Belgium) BVBA  
Honzebroekstraat 137 • 8800 Roeselare  
Phone: +32 +51 200 425  
Fax: +32 +51 200 449  
e-mail: [belgium@ultimed.org](mailto:belgium@ultimed.org)



November 2014-AL  
IA NB