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ПОЛОЖИТЕЛЕН:* Two distinct red lines appear. One line should be in the control region (C) and another line should be in the test region (T). A positive result indicates that H. pylori was detected in the sample.

* **NOTE:** The intensity of the red colour in the test line region (T) will vary depending on the concentration of H. pylori present in the sample. Therefore, any shade of red in the test region (T) should be considered positive.

NEGATIVE: One red line appears in the control region (C). No apparent red or pink line appears in the test region (T). A negative result indicates that H. pylori antigen is not present in the sample, or is present below the detectable level of the test.

INVALID: Control line fails to appear. Insufficient sample volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact distributor / manufacturer.

QUALITY CONTROL

A pink line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. If available, external controls may also be used to assure that the reagents are working properly and that the assay procedure is followed correctly. It is recommended that a control be tested at regular intervals as good laboratory testing process.

Users should follow the appropriate local guidelines concerning the running of external quality controls.

LIMITATIONS OF THE TEST

- The ulti med H. pylori Antigen Test Cassette is for professional in vitro diagnostic use, and should be used for the qualitative detection of Helicobacter pylori only.
- Following certain antibiotic treatments, the concentration of H. pylori antigens may decrease to the concentration below the minimum detection level of the test. Therefore, diagnosis should be made with caution during antibiotic treatment.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The ulti med H. pylori Antigen Test Cassette has been evaluated with specimens obtained from a population of symptomatic and asymptomatic individuals who presented for endoscopic examination.

Biopsy (Culture) served as the reference method for the H. pylori Antigen Test. Histology and a Rapid Urease Test (RUT) were performed on all negative culture specimens. The specimen was considered positive if Culture was positive. The specimen was also considered positive if the Culture was negative, but both Histology and RUT were positive.

H. pylori Antigen Test	Endoscope-based methods	
	+	-
+	132	0
-	0	154

Relative Sensitivity: >99,9% (97,3%~100,0%)*

Relative Specificity: >99,9% (97,6%~100,0%)*

Overall Agreement: >99,9% (98,7%~100,0%)*

*95% Confidence Interval

Analytical Sensitivity

The sensitivity of the ulti med H. pylori Antigen Test Cassette was determined by concentration of Urease. The sensitivity is 32 ng/ml of urease.

Specificity

Cross reactivity with following organisms has been studied at 1.0 x10⁹ organisms/ml. The following organisms were found negative when tested with the ulti med H. pylori Antigen Test Cassette:

Staphylococcus aureus	Proteus mirabilis	Neisseria gonorrhoea
Pseudomonas aeruginosa	Acinetobacter spp	Group B Streptococcus
Enterococcus faecalis	Salmonella choleraesuis	Proteus vulgaris
Group C Streptococcus	Gardnerella vaginalis	Enterococcus faecium
Klebsiella pneumoniae	Acinetobacter calcoaceticus	Hemophilus influenzae
Branhamella catarrhalis	E.coli	Neisseria meningitidis
Candida albicans	Chlamydia trachomatis	Rotavirus

REFERENCES

- Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
- Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
- Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
- Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:35S-41S.
- Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

Manufacturer	Contents sufficient for <n> tests
For in vitro diagnostic use only	Lot. no.
For single use only	Expiration date
Read instructions for use	Store at

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ulti med Products (Deutschland) GmbH